

Midas Rex® Onsite Planned Maintenance

Proactive + Preventive = Performance



Equipment availability and performance are key factors for an operating room. With Medtronic's Onsite PM service, a Medtronic factory-trained technician will come to your facility at a pre-arranged date and time. They will inspect and test your Midas Rex motors, attachments, consoles, and foot controls.

- Compliance: Your equipment is tested by a Medtronic factory-trained technician
- Convenience: Scheduled maintenance visit based on your schedule
- Availability: Onsite real-time exchanges may be performed*
- Dependability: Recommendations provided to your staff to help extend the life of your equipment.

Recommended Planned Maintenance Schedule

Maintenance Level	# of Tools Used per Motor per Year	Usage Rate per Motor per Year	PM Level
High	150-200	> 12	Every 3 months
Med High	100-149	8-11	Every 4 months
Med Low	50-99	4-7	Every 5 months
Low	< 50	< 4	Every 6 months

What's Included in a PM Visit

Inspection and testing is performed on each motor, console, foot-control, attachment, hose and cable by Medtronic factory-trained technician at your location.

- Visual, audible and tactile inspection for excessive wear and tear.
- Functional tests on motors, foot-controls and attachments with a factory supplied IPC.
- Functional tests on IPCs with factory supplied motors, foot-controls and attachments.
- Testing for several minutes with thermal and optical meters to validate proper operating specifications.
- Installation of latest software onto IPCs if applicable.
- A full PM inspection report for compliance record keeping.
- Recommendations for staff to ensure proper cleaning, handling, and sterilization.
- Equipment not passing inspection may be exchanged onsite or for next day delivery.



*Customers without an ESP will incur an additional charge

For more information contact your local Medtronic Neurosurgery area sales manager or visit www.responsecare.medtronic.com.

Medtronic plc
Neurosurgery
 4620 North Beach Street
 Fort Worth, Texas 76137
 USA
 Toll free: (800) 433-7080
 Telephone (817) 788-6400
www.medtronicneurosurgery.com



Caution: Federal law (USA) restricts this device for sale by or on the order of a physician. Refer to product instruction manual/package insert for instructions, warnings, precautions, and contraindications. Midas Rex® and ResponseCare® are registered trademarks of Medtronic, Inc.